

Precision – more than skin deep

SFM Spinal Needle Programme with Optimum Functionality



Spinal Needle with introducer

Different versions of hubs:



Needle hub SFM I



Needle hub SFM II

Tip geometry:



Pencil Point



Quincke

Product features

- Pencil Point Spinal Needles: Lateral opening located particularly close to the needle tip
- Quincke-type Spinal Needles: ("cutting" short bevel)
- Thin-walled tubes with large internal diameter enable advantageous cerebrospinal fluid (CSF) reflux
- Spinal Needles are also available in combination with introducer.
- Colour-coded stylet hub for easy and fast identification of size
- Highly transparent hub for easy control of reflux
- Delivery in customized packaging (EO sterilized) or in bulk

Gauge	External diameter Needle	Blade length
19	1,10 mm	70 mm - 90 mm -110 mm
20	0,90 mm	70 mm - 90 mm -110 mm
21	0,80 mm	70 mm - 90 mm -110 mm
22	0,70 mm	70 mm - 90 mm -110 mm
23	0,65 mm	70 mm - 90 mm -110 mm
24	0,55 mm	70 mm - 90 mm -110 mm
25	0,50 mm	70 mm - 90 mm -110 mm
26	0,45 mm	70 mm - 90 mm -110 mm
27	0,40 mm	70 mm - 90 mm -110 mm

Special sizes and minimum order quantities available on request

Please contact us:

Süddeutsche
Feinmechanik GmbH
Brückenstraße 5
D-63607 Wächtersbach

T +49 (0) 60 53/805-0
F +49 (0) 60 53/40 31
info@sfm.de
www.sfm.de

SFM Business Activities

Areas of Application and Products

- Anaesthesiology
(e.g. Spinal Needles, Epidural Needles, Epidural Sets)
- Biopsy
(Punches, Curettes)
- Blood Treatment
(e.g. Plasmapheresis Needles)
- Infusion Systems
(e.g. Strauss Needles, Port Needles)
- Neurology
(Thermolesion Needles)
- Cardiology
(Seldinger Needles)
- Transfer Systems
(e.g. Transfer Devices, Bottle Needles)
- Urology
(Split Cannulae)
- Veterinary Medicine
(Chip Cannulae, Implanter)

Services and Engineering

- Customer-specific development
- Tailor-made packaging solutions
- EO sterilisation in our own sterilisation facility
(validated according to EN 550:1994 and ISO 11135:1994)
- Certificates:
EN ISO 13485:2003
Permission to manufacture according to §13
(German Pharmaceutical Act)



Please contact us:

Süddeutsche
Feinmechanik GmbH
Brückenstraße 5
D-63607 Wächtersbach

T +49 (60 53) 805-0
F +49 (60 53) 40 31
info@sfm.de
www.sfm.de


medical devices